


|   |  |                           |             |
|---|--|---------------------------|-------------|
| <br>UNIVERSITY OF<br><b>TORONTO</b><br>MISSISSAUGA | <b>Chemical and<br/>Physical<br/>Sciences<br/>Department</b> | SOP #                     | <b>005</b>  |
|   |  | Revision #                | 02          |
|   |  | Implementation Date       | 2019-09-12  |
| Page #  | 1 of 4   | Last Reviewed/Update Date | 2019-09-12  |
| SOP Owner   | Grace Flock  | Approval                  | Grace Flock |

## Standard Operating Procedure: Liquid Nitrogen (N<sub>2</sub>) Handling

**WARNING:** LN<sub>2</sub> expands 700 times its volume displacing O<sub>2</sub> and creating and may cause asphyxiation. **IT DOES NOT SUPPORT LIFE**

**1. Purpose: Safe handling of Liquid Nitrogen**

---

**2. Scope:** Applies to all users

---

**3. Prerequisites: IRM training and RMF-SST**

---

**4. Responsibilities:** it is everybody responsibility to follow SOP and to report any equipment misuse and or deficiency to the research facility manager.

---

**5. Personal Protection Equipment (PPE): Closed toe shoes, Long Pants, Goggles, Faceshield**



**6. Procedure:**

---

### 6.1 Liq N<sub>2</sub> Storage:

Large amounts (more than 2 L) of LN<sub>2</sub> should always be stored in a well-ventilated area equipped with an O<sub>2</sub> sensor.



## 6.2 Handling/Dispensing LN<sub>2</sub>:

Transferring from primary container (Dewar) to large secondary containers (ex. 30L)  
These containers should be stored in the cryogenic room

- Put on PPE (as above pictograms)
- Always use the specially designed containers when transporting and handling LN<sub>2</sub> (see below for examples)



- Open valves of primary Dewar slowly to minimize thermal effects and control gas escape
- Do not fill secondary containers to more than 80% of capacity; expansion of gases may cause pressure buildup
- If the container tips over, let it go and evacuate (yourself and all people nearby).
- Following, contact your PI.

Bench top containers: utilized for small scale use in the labs/pods



- Never dispense liquid into an unapproved container, such as a Thermos® bottle. It will shatter.
- Transfer of LN<sub>2</sub> can cause splashing, wear long pants, face shield and appropriated gloves (provided on location, cryogenic storage room)
- Utilize specialized withdrawal devices instead of pouring (LN<sub>2</sub> Pump)
- Transfer liquid slowly to prevent thermal shock, pressure buildup, and splashing.
- Always wear appropriate PPE



### 6.3 Handling/Transporting/Thawing Cryotubes



Cryotubes used to contain samples stored under liquid nitrogen may explode without warning when handling and thawing.

When thawing cryotubes, take the following protective steps:

- Wear a face shield and safety goggles whenever handling cryogenic liquid.
- Wear appropriate insulated gloves
- Wear a buttoned lab coat and pants and closed toes shoes
- Place the cryotube in a secondary tube (example, falcon tube), as a shield, while transporting and or thawing

### 7.1 Emergency Response Procedure:

- 1) If there is a large spill or rupture of a Liq N<sub>2</sub> container, Evacuate. Spill may induce oxygen deficiency
- 2) Call 911 and Campus Police 905-569-4333.

## **Personal injuries:**

### Minor injuries:

- 1) Cold burns should be immediately flushed with tepid water or placed in a warm water bath.
- 2) Notify your PI.
- 3) DO NOT RUB SKIN – may damage tissue
- 4) Seek medical attention immediately for assessment and follow-up:  
**ex: emergency department at Credit Valley Hospital or HCC at UTM DV 1123A (depending on severity)**

### Major injuries: call 911

## **8. References:**

Prepared by Grace Flock (RMF manager)

Revised for University of Toronto Mississauga by Rocsana Pancescu (Chemistry Laboratory Coordinator, Davis Building, UTM Campus)

---